Lovibond® Water Testing

Tintometer® Group

Floc-Tester

ET830 • ET840 • ET850

(EN) Instruction manual..... 3	(DE) Gebrauchsanleitung..... 18
(ES) Manual de instrucciones.....33	(FR) Mode d'emploi..... 49
(PT) Manual de instruções.....65	(IT) Manuale d'istruzione..... 80
(NL) Gebruiksaanwijzing.....96	(TR) Kullanım Talimatları..... 111
(ZH) 使用说明..... 126	

Table of Contents

1 Introduction	4
1.1 General Information	4
1.1.1 Read instructions before use	4
1.1.2 Purpose of the Product	4
1.1.3 Authorized use	4
1.1.4 Requirements for safe use	4
1.1.5 Keep instructions	4
1.1.6 User qualification	5
1.1.7 Handling of hazardous chemicals	5
1.1.8 Disposal notes	5
1.2 List of all used signs in the document	5
2 Product overview	6
2.1 Guide to symbols	6
2.2 Certification	6
2.3 Features	7
2.4 Product description	7
2.4.1 Instrument view	7
2.4.2 Description of the switches	8
3 Commissioning	9
3.1 Operating environment	9
3.2 Contents of delivery	9
3.3 Installation	9
4 Operation	10
4.1 Control elements	10
4.1.1 List of control elements and their function	10
4.1.2 Buttons and function	10
4.2 Description of the operation procedures and the individual sections of the user interface	10
4.2.1 Overview	10
4.2.2 Mode MD1	11
4.2.3 Mode MD 2	11
5 Maintenance	13
5.1 Cleaning and maintenance	13
6 Troubleshooting	14
6.1 Errors and countermeasures	14
7 Accessories & Replacement parts	15
7.1 Spare parts list	15
8 Specifications	16

1 Introduction

1.1 General Information

1.1.1 Read instructions before use

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

1.1.2 Purpose of the Product

Lovibond Floc Testers provide you with optimum support when conducting your flocculation tests. They enable you to investigate the effect of different flocculants in different dosages under various test conditions on your sample material in order to determine the optimum parameters for flocculation in your process. Thanks to the various setting options, such as 2-stage programmes and variable lighting position (ET840 and ET850 only), every flocculation process can be explored in detail in a reproducible manner.

1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

1.1.4 Requirements for safe use

Note the following points for safe use:

- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- The product may only be used according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be used under the environmental conditions mentioned in this operating manual.
- The product must not be opened or modified.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storage conditions, see chapter "Specifications")

1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products. The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

1.1.8 Disposal notes

Dispose of the batteries and electrical devices at a suitable facility according to local legal requirements.

It is illegal to dispose of the batteries with household waste.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life.

Instruments marked with this symbol must not be disposed of in normal domestic waste.

1.2 List of all used signs in the document

The following symbols are used in this manual to identify sections that require special attention:

Danger!

A hazard exists that will result in death or severe injury if not avoided.

Notice!

Important information or specific instructions need to be strictly followed.

2 Product overview

2.1 Guide to symbols

Labels attached to the product should be strictly observed to avoid personal injury or damage to the product. Refer to this chapter for information regarding the nature of the danger or risk before taking any action where such label is present.

Do not use the device in rooms with a risk of explosion or in corrosive environments.

EN

Although the device offers a high level of protection against spilled liquids, it is necessary to disconnect the device from the power supply and clean it as a safety measure in the event of spilled liquids and splashes.

Ensure that the light screen is securely fastened, as it can be moved to two positions and may move during transport or when the device is tilted (ET840 and ET850 only).

This is a stirrer that moves masses. It is therefore important that the surface on which the device is placed is as flat and non-slip as possible so that the device can comply with EN IEC 61010-2-051.

The device has a safety coupling system to prevent damage to the mechanical drive if the stirrers are prevented from rotating during operation.

For professional users in the European Union:

If you wish to discard electrical and electronic equipment (EEE), please contact your dealer or supplier for further information.

For disposal in countries outside of the European Union:

This symbol is only valid in the European Union (EU). If you wish to discard this product please contact your local authorities or dealer and ask for the correct method of disposal.

2.2 Certification

EMC Immunity	EN 61326-1	CE mark
Safety	EN 61010-1	CE mark
	EN 61010-2-051	particular requirements
Further directives	2006/95/EC	low voltage equipment directive
	2004/108/EC	electromagnetic compatibility directive
	2006/42/EC	machinery directive

2.3 Features


Lovibond Floc Testers offer a high degree of user-friendliness. The many different setting options make it possible to determine the optimal conditions for flocculation processes in liquids.

- Four (ET830 & ET840) or six (ET850) stirring positions with individually adjustable stirrer height
- Suitable for beakers from 500 mL to 1500 mL (1500 mL on four stirring positions for the ET830 and on two or three stirring positions for the ET840 and ET850)
- Durability thanks to corrosion-resistant stainless steel stirrers
- Highly legible, illuminated display
- Easy-to-understand, language-independent operating concept
- Speed digitally adjustable in the range from 20 to 250 rpm
- Speed control ensures exact adherence to the selected test conditions
- Timer adjustable between 1 and 240 minutes and ∞ .
- Two independent speeds and timers with automatic changeover enable reproducible, 2-stage flocculation tests (e.g. short-term high speed for mixing, followed by slow stirring during coagulation over a longer period of time)
- Supports ASTM D-2035 compliant flocculation testing
- Robust handles for easy transport
- Portable operation with 12V DC possible (ET830)
- Lighting of the sample from behind or from below possible (ET840 and ET850)

EN

2.4 Product description

2.4.1 Instrument view



- 1 – Stirrer knob
- 2 – Display
- 3 – Keypad
- 4 – Stirrer paddle
- 5 – Power switch
- 6 – Lighting switch
- 7 – Transport handles
- 8 – Sample lighting

- 9 – Power cable connection
- 10 – Sample illumination fixing (ET840/ET850)
- 11 – Sample illumination connection (ET840/ET850)
- 12 – 12V DC connection (ET830)
- 13 – Fuse (ET830)

2.4.2 Description of the switches

EN

Switch	Function
	<p>Power switch Turns the device on and off</p>
	<p>Illumination switch Turns the sample illumination on and off</p>

3 Commissioning

3.1 Operating environment

The instrument can be used in any indoor or outdoor environment that is dry and deemed safe for an operator to perform the flocculation tests. This environment includes a temperature range from 5 to 50 degrees Celsius with a non-condensing humidity up to 80%. The instrument should be placed on a horizontal surface to achieve the best test performance.

3.2 Contents of delivery

- Floc tester
- Cable set
- Power adapter
- Car adapter (ET830)
- Operating instructions

3.3 Installation

Before connecting to the mains supply, check that the supply voltage matches the specification on the mains adapter. Failure to do so may result in unsafe operation, malfunction or damage to the device.

4 Operation

4.1 Control elements

4.1.1 List of control elements and their function

1	2	3	4
rprn	mode	timer	
v1: 2@	MD1	@m	
v2: 25@	off	@m	

5
6

1 – Stirring level v1
2 – Speed per stirring level
3 – Active mode
4 – Timer per stirring level
5 – Stirring level v2
6 – Stirring indicator

4.1.2 Buttons and function

Buttons	Functions
	Up button Increases the value at the cursor position
	Down button Decreases the value at the cursor position
	Start/Stop button Starts/stops the stirrer or the set stirring sequence
	Right/Mode button Press briefly to move the cursor one position further or to switch between the two stirring levels in MD1 mode Press and hold (3 seconds) to switch between MD1 and MD2 modes

4.2 Description of the operation procedures and the individual sections of the user interface

4.2.1 Overview

The device has two differently definable stirring levels, v1 and v2, and two different operating modes, MD1 and MD2. A speed can be preset for each stirring level. As soon as a stirring level is active, the stirrers rotate at the set speed. One timer is available for operating mode MD1 and two for operating mode MD2.

4.2.2 Mode MD1

In mode MD1, a stirring duration can be set using the timer for stirring level v2. This setting can be between 1 minute and 240 minutes. In addition, stirring for an unlimited period of time can be activated by setting this timer to 0. In mode MD1, the timer for stirring level just standing v1 displays the remaining stirring duration after stirring has started.

Pressing the start/stop button starts the preset mode. When the stirrers are started, the currently active values of a mode are marked with → arrows. Pressing the right/mode button in mode MD1 allows you to switch between the two preset stirring levels v1 and v2. The rotation speed of the currently active stirring level can be increased or decreased directly using the up and down buttons. The timer continues to run unchanged.

Figure 1:

	rpm	mode	timer
v1:	20←	MD1	->1m
v2:	240	on	1m

In Figure 1, the stirrers are running in MD1 mode at a speed of 20 rpm. The remaining stirring time is 1 minute. As a second stirring stage, 240 rpm can be selected here by briefly pressing the right/mode button (see Figure 2). An acoustic signal sounds when the timer expires.

Figure 2:

	rpm	mode	timer
v1:	20	MD1	->1m
v2:	240←	on	1m

In Figure 3, the stirrers are just standing. When the stirring programme starts, it is possible to switch between stirring speeds of 20 rpm and 250 rpm. The stirring duration is here set to unlimited stirring time.

Figure 3:

	rpm	mode	timer
v1:	20	MD1	@m
v2:	250	off	@m

4.2.3 Mode MD 2

In MD2 mode, both stirring levels run automatically one after the other. Each stirring level has its own timer.

Figure 4:

	rpm	mode	timer
v1:	120	MD2	1m
v2:	40	off	20m

In Figure 4, a speed of 120 rpm and a stirring time of 1 minute are set for stirring level v1, and a speed of 40 rpm and a stirring time of 20 minutes are set for stirring level v2. The stirrers don't run. By briefly pressing the right/mode button, the cursor can be moved to the timer for stirring level v1 and then set to 3 minutes, for example. The result can be seen in Figure 5:

Figure 5:

	rpm	mode	timer
v1:	120	MD2	3m
v2:	40	off	20m

Pressing the start/stop button starts the entire sequence of both stirring levels. When the first timer expires, an acoustic signal sounds and the second stirring level starts automatically.

Figure 6:

	rpm	mode	timer
v1:	120	MD2	0m
v2:	40	on	20m

The display then looks as shown in Figure 6. Note that the cursor has automatically jumped to the speed of stirring level v2. If necessary, this can be adjusted during the running sequence by pressing the up and down buttons. At the end of the sequence, an acoustic signal sounds again.

5 Maintenance

5.1 Cleaning and maintenance

To ensure that the height adjustment of the individual stirrers remains smooth, it is recommended to treat the shafts with Vaseline (never with grease or sprays) from time to time. Disconnect the device from the mains before starting cleaning. Clean at room temperature using water or isopropanol. Do not use any other organic solvents or aggressive cleaning agents, as these may damage the device.

6 Troubleshooting

6.1 Errors and countermeasures

Error description	Countermeasure
Device cannot be switched on	Check the power supply and fuse. When replacing fuses, make sure that only fuses with the same characteristics as the originals are used.
Individual stirrers do not run	The transmission belts may be worn (spare part 419183). Please contact our service department.
Stirrers sink down automatically	The stirrer holder is too loose or shows signs of wear. Please contact our service department.

EN

7 Accessories & Replacement parts

7.1 Spare parts list

Description	Item number
Operating instructions ET830/ET840/ET850	419175
Floc tester Drive set	419182
O-ring (ET830/ET840/ET850)	419183
Drive motor ET830/ET840/ET850	419814
LED strip ET830/ET840/ET850	419185
Main circuit board ET830/ET840/ET850	419186
Microchip with firmware ET830/ET840/ET850	419187

When ordering spare parts, please always state the serial number of your device in addition to the name of the device and the desired spare part. You will find this on the label on the back of the device.

EN

8 Specifications

Voltage		100 - 240 V
Frequency		50/60 Hz
Output Voltage Power Supply		15 VDC
Fuse		5x20 mm, F3.15AL250 V
Max. Power		50 W
Display		LCD backlit
Operation		Membrane keys
Power Supply		100 - 240 V / 50-60 Hz
Portability	ET830	Portable
	ET840	-
	ET850	-
Timer		0 - 240 minutes and ∞
Motor		DC Gear Motor with Encoder
Stirring Places	ET830	4 (portable)
	ET840	4
	ET850	6
Stirring Speed Control		20 - 250 revolutions per minute
Resolution		5 rpm
Precision		± 1 rpm
Rod and Stirrer Material		Stainless Steel
Environmental Conditions		5-50 °C
Stock Conditions		80 % max.rel. Humidity
IP Rating		IP 34
Compliance		CE
Protection Class		Power Supply: Class II, Device: Class III
Dimensions	ET830	255 x 428 x 253 mm
	ET840	800 x 410 x 200 mm
	ET850	
Weight	ET830	6.9 kg
	ET840	15.9 kg
	ET850	

Notice!

Subject to technical modification.

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Trademarks of the Tintometer Group
Tintometer Inc.